

HAEMOSTASIS



Company GmbH



INJECTION NEEDLES

HEMOCLIPS



INJECTION NEEDLES



ø 25 gauge

ø 23 gauge

ø 22 gauge

ø 19 gauge



Application of medications and solutions

- ✓ Automatic needle retraction
- ✓ Needle locking
- ✓ Colour coding




Company GmbH

mdd

Injection needles

Application of medications and solutions



Features

- ☐ Reliable needle guidance thanks to automatic needle retraction
- ☐ Rapid selection of the appropriate needle diameter thanks to permanent colour code
- ☐ Smooth puncturing of the target tissue
- ☐ Needle locking with extended or retracted needle
- ☐ Integrated ergonomic handle for comfortable use
- ☐ Individually sterile packed

| Needle diameter | Needle length | Outer catheter diameter | Working length | Pieces per box | Manufacturer Product Code |
|------------------------|---------------|-------------------------|----------------|----------------|---------------------------|
| FOR COLONOSCOPY | | | | | |
| 19 gauge / 1.0 mm | 4 mm | 2.4 mm | 230 cm | 10 | IN-23-19-4-230 |
| 19 gauge / 1.0 mm | 5 mm | 2.4 mm | 230 cm | 10 | IN-23-19-5-230 |
| 19 gauge / 1.0 mm | 6 mm | 2.4 mm | 230 cm | 10 | IN-23-19-6-230 |
| 22 gauge / 0.7 mm | 4 mm | 2.4 mm | 230 cm | 10 | IN-23-22-4-230 |
| 22 gauge / 0.7 mm | 5 mm | 2.4 mm | 230 cm | 10 | IN-23-22-5-230 |
| 22 gauge / 0.7 mm | 6 mm | 2.4 mm | 230 cm | 10 | IN-23-22-6-230 |
| 23 gauge / 0.6 mm | 4 mm | 2.4 mm | 230 cm | 10 | IN-23-23-4-230 |
| 23 gauge / 0.6 mm | 5 mm | 2.4 mm | 230 cm | 10 | IN-23-23-5-230 |
| 23 gauge / 0.6 mm | 6 mm | 2.4 mm | 230 cm | 10 | IN-23-23-6-230 |
| 25 gauge / 0.5 mm | 4 mm | 2.4 mm | 230 cm | 10 | IN-23-25-4-230 |
| 25 gauge / 0.5 mm | 5 mm | 2.4 mm | 230 cm | 10 | IN-23-25-5-230 |
| 25 gauge / 0.5 mm | 6 mm | 2.4 mm | 230 cm | 10 | IN-23-25-6-230 |



HEMOCLIPS



Hemoclip for endoscopic marking and hemostasis

- ✓ Secure clip placement
- ✓ 360° rotation
- ✓ Colour coding

mdd

Heliclip

Hemoclip for endoscopic marking and hemostasis



Features

- ☐ 180 degree opening and stepless 360° rotation of the clip facilitate placement in the tissue
- ☐ Safe and easy release of the clip
- ☐ PTFE tube coating ensures optimum sliding properties
- ☐ Atraumatic tip protects the endoscope from damage
- ☐ Ergonomic handle supports a comfortable application
- ☐ Colour-coded for quick orientation in choosing the appropriate clip length
- ☐ Individually sterile packed
- ☐ Available from March*, June**, July 2020***

| Clip length | Opening | Outer catheter diameter | Working length | Pieces per box | Manufacturer Product Code |
|-------------|---------|-------------------------|----------------|----------------|---------------------------|
| 9 mm*** | 18 mm | 1.8 mm | 230 cm | 10 | HC-23-230-9 |
| 12 mm** | 24 mm | 1.8 mm | 230 cm | 10 | HC-23-230-12 |
| 16 mm* | 32 mm | 1.8 mm | 230 cm | 10 | HC-23-230-16 |



Company GmbH



UpViser Aps
Adelgade 45, 1
8660 Skanderborg, Denmark
Tel.: +45 6171 7625
E-Mail: salg@upviser.com
www.upviser.dk



Exclusive Partner Europe
Body Products GmbH

mdd Company GmbH
Rengglochstr. 38
CH-6012 Obernau, Switzerland
Tel: +39 348 0026970
E-Mail: gb@mddcompany-gmbh.com
www.mddcompany-gmbh.com