



## PRODUCT SAFETY DATA SHEET

Salts Healthcare Limited  
Richard Street Aston Birmingham B7 4AA  
Telephone: +44 (0) 121 333 2000  
Fax: +44 (0) 121 333 2010  
E-mail: salt@salts.co.uk  
Website: www.salts.co.uk

**PRODUCT NAME(S):**

WipeAway Adhesive Remover

**2. INTENDED USES:**

For adhesive removal and skin cleansing. For EXTERNAL USE ONLY

**3. COMPOSITION:**

A flexible laminate sachet containing a disposable non-woven swab impregnated with a non aqueous solution of Iso-Alkane Mixture, Emollient Oil and fragrance. ALCOHOL FREE

**4. PHYSICAL AND CHEMICAL PROPERTIES:**

WipeAway wipes are flammable (flashpoint  $>56^{\circ}\text{C}$ ).

**5. HEALTH HAZARD:**

No health hazard is anticipated during normal handling of this product. Avoid contact with eyes - may cause irritation.

**6. FIRE HAZARD AND EMERGENCY ACTION:**

Flammable. Packaging is a combustible solid which gives off toxic fumes, mainly carbon monoxide and carbon dioxide, when ignited. Extinguish fires with water fog, foam, dry chemicals, halons or carbon halons or carbon dioxide. Wear self-contained breathing apparatus.

**7. STORAGE PRECAUTIONS:**

Store in a cool dark place. Keep away from sources of ignition - No Smoking.



**PRODUCT SAFETY DATA SHEET –  
WIPEAWAY ADHESIVE REMOVER  
Page 2 of 2**

**SALTS**  
HEALTHCARE

**8. TRANSPORT PRECAUTIONS:**

Product is classed as flammable.

**9. HANDLING/USE/PROTECTIVE CLOTHING:**

Withdraw from sachets only the correct number of wipes required.

**10. DISPOSAL:**

Discard used wipes safely away from sources of ignition. Controlled disposal according to local environmental health, flammable hazardous waste disposal guidelines.

**11. FIRST AID:**

- a) Inhalation: Large quantities may cause coughing or wheezing.
- b) Contact with skin: Low irritation in normal individuals.
- c) Contact with eyes: May cause irritation, flush with copious quantities of water.
- d) Ingestion May cause lung damage if large quantities are swallowed.

**12. ADDITIONAL INFORMATION:**

Avoid strong oxidising agents

**13. NAME, ADDRESS AND TELEPHONE NUMBER OF SUPPLIER**

Salt & Son Ltd  
Lord Street  
Birmingham  
B7 4DS

TEL: 0121 333 2000  
FAX: 0121 359 0830

**14. REFERENCE NUMBER, DATE OF ISSUE:**

WA1

NOVEMBER 2007

This information is provided in accordance with the requirements of the UK Health and Safety at Work Act 1974 and specifically in order to assist risks' of the product to make their 'Assessment of health risks' as required by the UK Control of Substances Hazardous to Health Regulation 1988 (COSHH assessments). Provision of this information does not preclude users from seeking advice from other sources as indicated in the COSHH guides.

This information is intended to cover potential hazards at the place of work and does not detail medical uses, indications, contra-indications and precautions for the treatment of patients.  
Q:\COSHH\WipeAway I2.doc

