Corsinel®

BELT W/PANEL


Prepared for individual hole




Unique panel with option for individual hole


4-way stretch for controlled compression


Optimal comfort and freedom to move



Art. no. 3412







Prepared for cutting of hole template included

The Corsinel abdominal belt with panel is designed with 4-way stretch offering controlled compression. This helps prevent development of hernia.

The uniquely designed front panel offers the opportunity to cut a hole to fit the individual stoma pouch, while containing the controlled compression of the stoma/hernia.

The seamless belt conforms to the body contour, providing optimal comfort and freedom to move.

Indications:

Medical device for support and compression for stoma, parastomal hernia, hernia treatment or prevention, and where an individual hole for the stoma pouch may be needed.

Precautions:

Always consult a doctor or a nurse with questions regarding use or correct fitting. Make sure that the stoma pouch is placed correctly underneath the belt in order to avoid pressure points and/or obstructed flow from the stoma. If the individual hole is not cut to size in accordance with the instructions, or the belt is incorrectly used, there is the risk that occlusion, entrapment, prolapse, leakage or discomfort may occur.

If any unexpected discomfort, stoma or skin irritation or similar occur, please consult your doctor or nurse immediately.

Height of belt* (front)

The belt may feel warm to use as it is very tight-fitting, however, this is no cause for concern, but consult your doctor or nurse in case of discomfort or skin irritation. No instances of discomfort caused by heat or moisture build-up have been reported or noticed by patients or health care professionals.

Size chart


3XI

SIZE СМ INCH СМ INCH 78 - 90 31 - 35 21,5 81/2 М 86 - 98 34 - 3922.5 9 25.5 10 94 - 106 37 - 42102 - 116 40 - 46 28,0 11 44 - 51 112 - 130 28,0

50 - 58

126 - 148 Front/unstretched +/- 1 cm / 1/2 inch

Recommended size chart (abdomen)



Composition:

50% polyamide 50% elastane

Washing instructions:



Close the Velcro in the outer position before washing. Wash with similar colours and in a washing bag.

We recommend wash at 60°C. The belt can be washed at 95°C, but it reduces the durability of the belt.

Washing or drying at temperatures higher than recommended will decrease durability and function. Use of chlorine-based bleach or chlorine-based disinfectant will decrease durability and function.

For further information please contact:







