
	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 1	Page 1 of 1
	Issue date: 04.04.2022	Effective date: 04.04.2022

SUCTION SET

REF. NO. **1128 3532 80 ML**



COMPONENT DESCRIPTION

COMPONENT NO.

1. SUCTION PIECE

rigid tube ID/OD 3.0/4.0, tube on funnel 186 mm, bent
material: PVC

0.28018640

2. HANDLE SHORT

handle white ID 4.0 mm, 160 mm length with **vacuum control**
material: Styrolux (SBC)

4.00016040

3. TUBE

Multi-Lumen transparent tube ID/OD 5.6/8.0 mm,
material: PCV 73°ShorA, DEHP – FREE,
• 1128 3532 80 ML – length 3400 mm

0.56340080

4. PROXIMAL END

funnel (funnel) 8.0 PVC
material: PCV 50°ShorA, DEHP – FREE, blue

0.00000080

PACKAGING MATERIAL

Innerbag (3500) 180 x 380 mm

2.18038004

Peelpack (packing machine) 200 x 350 mm

5.20035410

Outerbag 310 x 600 mm

25 pcs/bag

2.31060010

Box 386 x 386 x 260 mm

50 pcs/box

6.38638626

Label on box 85 x 131 mm

7.00085131

Work instruction

QI 01

Packaging instruction

QI 04

Comments:

	Glued	Not glued
Suction piece and handle	X	
Suction piece and handle short put into innerbag.		