
	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 1	Page 1 of 1
	Issue date: 04.04.2022	Effective date: 04.04.2022

SUCTION SET

REF. NO. 1148 2532 80 ML



COMPONENT DESCRIPTION

COMPONENT NO.

1. SUCTION PIECE

rigid tube ID/OD 3.0/4.0, 186 mm, 4 distal eyes, bent
material: PVC

0.48022040

2. HANDLE SHORT

handle white ID 6.5 mm, 160 mm length, with **vacuum control**
material: Styrolux (SBC)

4.20016065

3. TUBE

Multi-Lumen transparent tube ID/OD 5.6/8.0 mm,
material: PVC 73°ShoreA, DEHP – FREE,
• 1148 2532 80 ML – 3400 mm length

0.56240080

3. PROXIMAL END

funnel 8.0 PVC
material: PVC 50°ShoreA, DEHP – FREE, blue

0.00000080

PACKAGING MATERIAL

Innerbag 180 x 380 mm		2.18038004
Peelpack (packing machine) 200 x 350 mm		5.20035410
Outerbag 310 x 600 mm	25 pcs/bag	2.31060010
Box 386 x 386 x 260 mm	50 pcs/box	6.38638626
Label on box 85 x 131 mm		7.00085131

Work instruction	QI 01
Packaging instruction	QI 04

Comments:

	Glued	Not glued
Suction piece and handle	X	
Suction piece and handle short put into innerbag.		