



EYEMAX

SINGLE-USE VIDEO CHOLANGIOPANCREATICOSCOPY









HIGH-PERFORMANCE PROCESSOR

EFFICIENT IMAGE PROCESSING

- Qualitative transfer of video images with high resolution
- User-friendly operation
- Image enlargement up to 3x
- Adjustable display of the visualisation

EXCELLENT IMAGE QUALITY

The EYEMAX cholangiopancreaticoscope enables excellent high resolution image quality for improved visualisation of the bile and pancreatic ducts which can reduce the number of examinations.



MICRO-TECH 9

HARNESS

FOR MAXIMUM USER COMFORT

- Ergonomic alternative for conventional secure fixing
- Relief during longer examinations as the weight of the cholangiopancreatico- and duodenoscope do not have to be held in one hand
- One-hand operation of the cholangiopancreaticoscope guaranteed



Secure fixing is achieved by the anchor on the handle of the cholangiopancreaticoscope and attachment to the harness.



The cholangiopancreaticoscope is latched and secured to the harness to enable one-hand operation.



COMPATIBLE ACCESSORIES

FOR SUCCESSFUL OPERATION

- Large working channel with 1.2 mm or 2.0 mm diameter for conclusive samples and easier sampling
- Easy insertion of forceps even with a bent endoscope due to the short, rigid distal end of the forceps
- In addition to the irrigation channels, fluids can also be injected via the working channel



SPECIFICATIONS

Packing unit: 1 piece

REF	Max. catheter Ø mm	Wor	Working length mm		
CHOLANGIOPANCRE	ATICOSCOPE				
CDS11005	3.9		2.0	650	
CDS11006	3.9		2.0	2200	
CDS22005	3.2		1.2	650	
CDS22006	3.2		1.2	2200	
Packing unit: 1 piece					
REF	Interfaces	Power supply	Dimensions (LxWxH) mm	mm Weight	
PROCESSOR					
BS-W-100	DVI-D, Composite	230V	355 x 297 x 107	6 kg	
Packing unit: 1 piece					
REF					
HARNESS					
409501002					

REF	Min. working channel Ø mm	Length mm	Opening width mm	Jaw Ø mm	Spike	Covering	Туре
ACCESSORIES (FORCEPS)							
BF10001	1.1	1200	3	1.0	Without	PTFE	Alligator
BF10006	1.1	2900	3	1.0	Without	PTFE	Alligator
BF16001	1.8	1200	4.5	1.6	Without	PE	Alligator
BF16006	1.8	2900	4.5	1.6	Without	PE	Alligator
Packing unit: 2 pieces							

10 MICRO-TECH 11



