





Product data sheet - 15.45.10

24.08.2015

General Data			
Included in	<b>Catalogues</b> 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.60 (Plastic neu),		
Last update	2011-05-03		
Title	<b>Hemostatic Forceps</b>		
Article number	<b>15.45.10</b>	Author	Ultra
Description	Hemostatic Forceps, super delicate, highest precision, light weight. Ideal for: Vascular, Hand, Micro, Nuero and Plastic Surgery		
Length	10 cm - 4"	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826023133		
GMDN	34948		
Risk category EC	1		
UMDNS	10-865		
Net weight	0.0080		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 15.45.10

24.08.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---