





Product data sheet - 213.03.14

24.11.2015

General Data			
Included in	Catalogues 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.54 (OMF Instr), 450.00.58 (Cranialis), 450.00.59 (Spinalis), Brochures 451.03.23 (Medicut),		
Last update	2008-01-24		
Title	Operating Scissors		
Article number	213.03.14	Author	black line
Description	Operating Scissors		
Length	14.5 cm - 5 3/4"	Weight	
		Diameter	
Properties	curved	Tip / Working end	blunt/blunt
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826173494		
GMDN	38727		
Risk category EC	1		
UMDNS	13-480		
Net weight	0.0450		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 213.03.14

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material	T_GOLD_T_CHROM	Teil eines Instrumentes ist Vergoldet weiteers Teil verchromt zB Schraube Ringe	

Delivery mode	
---------------	--