




Product data sheet - 52.31.03

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu),		
Last update	2013-10-15		
Title	<b>Uterine Dilator</b>		
Article number	<b>52.31.03</b>	Author	Hegar
Description	Uterine Dilator		
Length		Weight	
		Diameter	Ø 3,5 mm
Properties	single-ended	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826066147		
GMDN	11267		
Risk category EC	1		
UMDNS	11-266		
Net weight	0.0150		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 52.31.03

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	-------------------------------------------------------------------------------------