





Product data sheet - 40.25.14

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.50 (Ortho), 450.00.63 (Ortho),		
Last update	2014-04-28		
Title	<b>Bandage Scissors</b>		
Article number	<b>40.25.14</b>	Author	Lister standard
Description	Bandage Scissors		
Length	14 cm - 5 1/2	Weight	
		Diameter	
Properties		Tip / Working end	angled to side
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826044404		
GMDN	13481		
Risk category EC	1		
UMDNS	13-481		
Net weight	0.0470		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 40.25.14

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---