





Product data sheet - 202.21.61

24.11.2015

General Data			
Included in	Catalogues 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.60 (Plastic neu), Brochures 451.03.23 (Medicut),		
Last update	2008-01-24		
Title	Iris Scissors		
Article number	202.21.61	Author	supercut plus
Description	Iris Scissors		
Length	11.5 cm - 4 1/2	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826171230		
GMDN	13488		
Risk category EC	1		
UMDNS	13-488		
Net weight	0.0180		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 202.21.61

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material	WSG_20-GO-55_CRTZ	RCoCr-C	
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode	
---------------	--