





Product data sheet - 26.35.02

24.11.2015

General Data			
Included in	--		
Last update	2014-02-03		
Title	Raspatory		
Article number	26.35.02	Author	Freer
Description	Raspatory		
Length	19 cm - 7 1/2"	Weight	
		Diameter	
Properties		Tip / Working end	5 x 5 mm, sharp/blunt, double-ended
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	<i>No preparation instructions available.</i>		
Operation manuals	<i>No operation manuals available.</i>		
GTIN			
GMDN	35048		
Risk category EC	1		
UMDNS	15-218		
Net weight	<i>No entries available.</i>		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 26.35.02

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---