





Product data sheet - 24.90.30

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.48 (Plastic), 450.00.57 (Plastic), 450.00.60 (Plastic neu),		
Last update	2011-02-25		
Title	<b>Areola Marker</b>		
Article number	<b>24.90.30</b>	Author	
Description	Areola Marker, blunt. For precise marking of the nipple during combined breast-lifting and reduction of the areola.		
Length		Weight	
		Diameter	Ø 42 mm
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826034641		
GMDN	32313		
Risk category EC	1		
UMDNS	12-443		
Net weight	0.0200		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 24.90.30

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Turned part material			
Handle material			
Others parts material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Coating material			

Delivery mode

