





Product data sheet - 57.28.53

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.37 (ENT neu), 450.00.58 (Cranialis), 450.00.59 (Spinalis), <b>Brochures</b> 451.04.32 (minimized), 451.04.56 (Kerrison),		
Last update	2014-08-08		
Title	<b>Dismantable Laminectomy Rongeurs "DieBlauen." with silicone coated handles</b>		
Article number	<b>57.28.53</b>	Author	magic cut Kerrison
Description	Dismantable Laminectomy Rongeurs "DieBlauen." with silicone coated handles. Magic Cut is a special coating method to obtain a considerably higher surface hardness.		
Length	18 cm - 7 1/8"	Weight	
		Diameter	
Properties	narrow footplate, spring golden, one piece dismantable rongeur, firmly integrated ejector, external guidance, jaw:13 mm	Tip / Working end	40° <b>up-cutting</b> , jaws 3 mm wide
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	G572800_2014_07		
GTIN	4046826404451		
GMDN	35287		
Risk category EC	1		
UMDNS	13-231		
Net weight	0.1710		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 57.28.53

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Handle material	ELASTOSIL LR3003		
Others parts material			
Coating material	TIALN	Instrument ist Titan Alu Nitrit beschichtet zB Magic Cut	

Delivery mode	
---------------	--