





Product data sheet - 13.80.07

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu),		
Last update	2013-09-20		
Title	Suture Clips		
Article number	13.80.07	Author	Michel
Description	Suture clips		
Length	7.5 x 1.75 mm	Weight	
		Diameter	
Properties	stainless steel	Tip / Working end	
Article link		Sales unit	Unit for sale = 1000 pieces (10 packages with 100 each)
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826020392		
GMDN	38143		
Risk category EC	2a		
UMDNS	15-964		
Net weight	0.0880		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 13.80.07

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---