





Product data sheet - 15.45.12

26.11.2015

General data			
Contained in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.54 (OMF Instr), 450.00.57 (Plastic), 450.00.60 (Plastic neu), <b>Brochures</b> 451.03.00 (Plastisch), 451.03.70 (Nagata),		
Last update	2011-05-03		
Title	<b>Hemostatic Forceps</b>		
Article number	<b>15.45.12</b>	Author	Halsted-Mosquito
Description	Hemostatic Forceps		
Length	12.5 cm - 5	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Numbers / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	• MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826023140		
GMDN	34948		
Risk category EC	1		
UMDNS	10-865		
Net weight	<i>No entries available.</i>		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 15.45.12

26.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode

