





Product data sheet - 214.32.14

24.11.2015

General Data			
Included in	Catalogues 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), Brochures 451.03.23 (Medicut),		
Last update	2008-01-24		
Title	Dissecting Scissors		
Article number	214.32.14	Author	Metzenbaum black line
Description	Dissecting Scissors		
Length	14 cm - 5 1/2"	Weight	
		Diameter	
Properties	straight	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826173838		
GMDN	13492		
Risk category EC	1		
UMDNS	13-492		
Net weight	0.0290		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 214.32.14

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material	T_GOLD_T_CHROM	Teil eines Instrumentes ist Vergoldet weilers Teil verchromt zB Schraube Ringe	

Delivery mode

