





Product data sheet - 214.33.20

24.11.2015

| General Data | | | |
|---|---|---|-----------------------|
| Included in | Catalogues 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), Brochures 451.03.23 (Medicut), | | |
| Last update | 2008-01-24 | | |
| Title | Dissecting Scissors | | |
| Article number | 214.33.20 | Author | Metzenbaum black line |
| Description | Dissecting Scissors | | |
| Length | 20 cm - 8" | Weight | |
| | | Diameter | |
| Properties | curved | Tip / Working end | |
| Article link | | Sales unit | |
| Accessories | | Number / Size | |
| Complete picture | | Detail images | |
|  | |  | |
| Quality management data | | | |
| Preparation instructions | MPK02 | | |
| Operation manuals | <i>No operation manuals available.</i> | | |
| GTIN | 4046826173951 | | |
| GMDN | 13492 | | |
| Risk category EC | 1 | | |
| UMDNS | 13-492 | | |
| Net weight | 0.0560 | | |



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 214.33.20

24.11.2015

| Materials | | | |
|-------------------------|-----------------|--|--------------------------------|
| | Material number | Description | Elements |
| Main component material | 1.4117 | X38CrMoV15 | C 0,35 - 0,40% Cr 14 - 15% |
| Turned part material | 1.4305 | X8CrNiS18-9 | Cr 17,0 - 19,0% Ni 8,0 - 10,0% |
| Handle material | | | |
| Others parts material | | | |
| Coating material | T_GOLD_T_CHROM | Teil eines Instrumentes ist Vergoldet weilers Teil verchromt zB Schraube Ringe | |

Delivery mode

