





Product data sheet - 04.13.13

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.54 (OMF Instr), 450.00.57 (Plastic), 450.00.60 (Plastic neu),		
Last update	2008-01-24		
Title	Dissecting Scissors		
Article number	04.13.13	Author	Ragnell (Kilner) standard
Description	Dissecting Scissors		
Length	13 cm - 5 1/8	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826003760		
GMDN	38727		
Risk category EC	1		
UMDNS	13-480		
Net weight	0.0290		



Product data sheet - 04.13.13

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	--