





Product data sheet - 29.85.14

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.50 (Ortho), 450.00.63 (Ortho),		
Last update	2009-03-19		
Title	Sharp Curette, oval		
Article number	29.85.14	Author	Volkman
Description	Sharp Curette, oval		
Length	14.5 cm - 5 3/4	Weight	
		Diameter	
Properties		Tip / Working end	double-ended
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826042431		
GMDN	31335		
Risk category EC	1		
UMDNS	11-084		
Net weight	0.0230		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 29.85.14

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---