





Product data sheet - 29.83.12

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.50 (Ortho), 450.00.63 (Ortho),		
Last update	2009-03-19		
Title	Sharp Curette, oval		
Article number	29.83.12	Author	Williger
Description	Sharp Curette, oval		
Length	14 cm - 5 1/2	Weight	
		Diameter	
Properties		Tip / Working end	double-ended
Article link		Sales unit	
Accessories		Number / Size	# 1-2
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826042332		
GMDN	31335		
Risk category EC	1		
UMDNS	11-084		
Net weight	0.0100		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 29.83.12

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---