





Product data sheet - 19.16.26

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.54 (OMF Instr),		
Last update	2008-01-25		
Title	Dressing Forceps		
Article number	19.16.26	Author	Maier
Description	Dressing Forceps		
Length	26 cm - 10 1/4	Weight	
		Diameter	
Properties	straight, with catch, box-lock	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826027230		
GMDN	11791		
Risk category EC	1		
UMDNS	11-791		
Net weight	0.1010		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 19.16.26

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---