





Product data sheet - 88.17.06

24.11.2015


| General Data | | | |
|---|--|---|--------|
| Included in | Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.45 (Neuro), 450.00.58 (Cranialis), 450.00.59 (Spinalis), | | |
| Last update | 2013-09-23 | | |
| Title | Coagulation Forceps, bipolar, straight, blunt | | |
| Article number | 88.17.06 | Author | |
| Description | Coagulation Forceps, bipolar, straight, blunt | | |
| Length | 16,5 cm - 6 1/2" | Weight | |
| | | Diameter | |
| Properties | sterelizable insulating plastic coating | Tip / Working end | 0.6 mm |
| Article link | | Sales unit | |
| Accessories | | Number / Size | |
| Complete picture | | Detail images | |
|  | |  | |
| Quality management data | | | |
| Preparation instructions | MPK02 | | |
| Operation manuals | G880702_2014_12 | | |
| GTIN | 4046826163297 | | |
| GMDN | 47848 | | |
| Risk category EC | 2b | | |
| UMDNS | 11-502 | | |
| Net weight | 0.0310 | | |



Product data sheet - 88.17.06

24.11.2015

| Materials | | | |
|-------------------------|-----------------|-------------|----------------------------|
| | Material number | Description | Elements |
| Main component material | 1.4021 | X20Cr13 | C 0,16 - 0,25% Cr 12 - 14% |
| Turned part material | | | |
| Handle material | | | |
| Others parts material | | | |
| Coating material | RILSAN | | |

| | |
|---------------|---|
| Delivery mode |  |
|---------------|---|