





Product data sheet - 11.70.22

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.53 (Cardio neu), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), Brochures 451.03.10 (Sch,Pinz,Nh),		
Last update	2014-01-30		
Title	Micro Needle Holder		
Article number	11.70.22	Author	Castroviejo-micro
Description	Micro Needle Holder		
Length	14.5 cm - 5 3/4	Weight	
		Diameter	
Properties	straight, without catch, smooth	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826018009		
GMDN	12726		
Risk category EC	1		
UMDNS	15-621		
Net weight	0.0560		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 11.70.22

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---