




Product data sheet - 47.20.27

24.11.2015

| General Data | | | |
|---|--|-------------------|--|
| Included in | Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), | | |
| Last update | 2011-10-27 | | |
| Title | Kidney Bowl | | |
| Article number | 47.20.27 | Author | |
| Description | Kidney Bowl | | |
| Length | | Weight | |
| | | Diameter | |
| Properties | 270 mm, 18/8 stainless steel | Tip / Working end | |
| Article link | | Sales unit | |
| Accessories | | Number / Size | |
| Complete picture | | Detail images | |
|  | | | |
| Quality management data | | | |
| Preparation instructions | MPK01 | | |
| Operation manuals | <i>No operation manuals available.</i> | | |
| GTIN | 4046826050689 | | |
| GMDN | 42893 | | |
| Risk category EC | 1 | | |
| UMDNS | 13-692 | | |
| Net weight | 0.3450 | | |




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 47.20.27

24.11.2015

| Materials | | | |
|-------------------------|-----------------|-------------|--------------------------------|
| | Material number | Description | Elements |
| Main component material | 1.4301 | X5CrNi18-10 | Cr 18,0 - 19,5% Ni 8,0 - 10,5% |
| Turned part material | | | |
| Handle material | | | |
| Others parts material | | | |
| Coating material | | | |

| | |
|---------------|---|
| Delivery mode |  |
|---------------|---|