



# Stand scale with height rod

7835

- ▶ Automatic BMI calculation
- ▶ Comfortable display height
- ▶ Robust platform
- ▶ Versatile functionalities
- ▶ Mobile use



## Features

The stand scale is the ideal scale for fitness studios, pharmacies, care or other social institutions. Due to the raised display, the weight value can be read off comfortably at eye level and through illuminated digits. The tread surface is robust and easy to clean. The scale is a special highlight due to its various functionalities. For example, with the help of the hold function, the weight value can also be read off only after weighing. The mother-child weighing function enables the weight of small children who cannot yet stand on the scales themselves to be determined. The scale also calculates the BMI index. After entering your height, you will receive your BMI index immediately during weighing. With the integrated interface the values can be transmitted directly to a PC or printer.

The scale can be operated by a power supply unit or be used with batteries for mobile applications. The automatic switch-off function switches the scale off when not in use, thus ensuring that no unnecessary energy is consumed. Integrated castors also make it easy to transport the scale.

- ▶ Including length measuring rod, unit: 1 mm
- ▶ Automatic calculation of the BMI by entering the body height and the determined weight value
- ▶ Display of weight, height and BMI
- ▶ Large tread area: 350 x 350 mm
- ▶ Brightly illuminated LC display with 20 mm high digits
- ▶ Easy to clean surface
- ▶ Taring function: documents placed on the machine can be tared
- ▶ Hold function for easy reading of the weight value after weighing
- ▶ Mother and child weighing
- ▶ Interface RS232 optional, for connection of a printer or PC
- ▶ Unit changeover kg/lb
- ▶ Energy-saving Auto-Off function
- ▶ Battery operation with 6 AA batteries (included)
- ▶ Power supply operation optional (see accessories)
- ▶ Integrated rollers for easy transport of the scale to the patient

## Products

Order number Variation	Approvable	Division g	Maximum load kg
7835.02.001 Stand scale with height rod	No	100	200
7835.02.002 Stand scale with height rod	No	100	300

## Accessories

Order number Variation	Description
2557.05.002	Power supply unit 100 - 240 V, 50-60 Hz
2557.25.003	UK exchange adapter
2550.03.002	Data cable (cable length 2 m)
2563.41.005	Interface RS 232, installed (must be included in the purchase order, no retrofit possible)