



# BASKET EXTRACTOR

## REMOVE FOREIGN OBJECTS SAFELY AND EFFECTIVELY

The extraction baskets from Bauer medical are not only suitable for stone extraction in the bile duct, but also for use in the upper and lower intestinal tract for the safe recovery

and removal of foreign bodies such as marbles, batteries or mouth fitting toys. The rotation of the Vario basket ensures safe and reliable recovery of the foreign body.

### SPECIFIC CHARACTERISTICS

- Fixed pre-determined breaking point
- Smooth wire with optimized gliding quality
- Ergonomic grip for improved tactile feeling
- Stone extractor basket with atraumatic tip
- Vario basket



Atraumatic tip gentle to working channel and mucosa



Highly mobile grip for precision insertion and extraction

## SPECIFICATIONS

REF	Form	Features	Ø basket mm
<b>DISPOSABLE BASKETS, PULL ROPE-LENGTH 200 CM</b>			
BK1-B10-25-23-200	Vario	Multifilament, 150°-rotation, protection head	15 / 25
BK1-B10-30-23-200	Vario	Multifilament, 150°-rotation, protection head	20 / 30
BK1-B20-20-23-200	Hexagonal	Multifilament, 4 rope, protection head	20
BK1-B20-25-23-200	Hexagonal	Multifilament, 4 rope, protection head	25
BK1-B20-30-23-200	Hexagonal	Multifilament, 4 rope, protection head	30
BK1-B20-35-23-200	Hexagonal	Multifilament, 4 rope, protection head	35
BK1-B21-35-23-200	Hexagonal	4 rope, half coiled	35
BK1-B22-20-23-200	Hexagonal	4 rope, full coiled	20
BK1-B22-25-23-200	Hexagonal	4 rope, full coiled	25
BK1-B22-30-23-200	Hexagonal	4 rope, full coiled	30
BK1-B22-35-23-200	Hexagonal	4 rope, full coiled	35
<b>DISPOSABLE GUIDE WIRE BASKETS, PULL ROPE-LENGTH 200 CM</b>			
BK1-B90-25-27-200	Combi-hexagonal	Multifilament, retractable tip	15 / 25
<b>DISPOSABLE LITHOTRIPSY BASKETS, PULL ROPE-LENGTH 400 CM</b>			
LT1-B10-25-23-400	Vario	Multifilament, 150°-rotation, retractable tip	15 / 25
LT1-B10-30-23-400	Vario	Multifilament, 150°-rotation, retractable tip	20 / 30
LT1-B20-20-23-400	Hexagonal	Multifilament, retractable tip	20
LT1-B20-25-23-400	Hexagonal	Multifilament, retractable tip	25
LT1-B20-30-23-400	Hexagonal	Multifilament, retractable tip	30
LT1-B20-35-23-400	Hexagonal	Multifilament, retractable tip	35

Packaging unit: 5 pieces