

Targeted Temperature Management

KOOL-KIT®

The Kool-Kit® includes a Head Wrap, Patient Vest and Lower Body Blanket to maximize coverage and temperature accuracy. The Kool-Kit® provides caregivers ready access to their patients, and helps administer temperature management and controlled re-warming for each patient. This method of cooling allows for early treatment that is fast, cost-effective and user-friendly.

Notice that the Head Wrap and Patient Vest were designed to cover the subclavian vein, subclavian artery and superior vena cava, all of which do not constrict with mild hypothermia.

- About 7% of blood flows through the area treated by the head wrap every minute.¹
- About 33% of blood flows through the vasculature exclusively treated by the Kool-Kit® vest.¹
- About 26% of blood flows through the groin area treated by the lower body blanket.¹
- Can cover approximately 66% of your patient, allowing for excellent body surface coverage.¹

1. Data on file at Gentherm



Kool-Kit®

TARGETED TEMPERATURE MANAGEMENT

Ordering Information

Cat No.	Description	Quantity
900	Kool-Kit® (1 head wrap, 1 patient vest & 1 lower body blanket)	5/box
910	Kool-Kit® Head Wrap (1 head wrap & 1 lower body blanket)	5/box
920	Kool-Kit® Patient Vest (1 patient vest & 1 lower body blanket)	5/box
930	Large Kool-Kit® (1 head wrap, 1 patient vest, 1 full body blanket)	5/box
950	Kool-Kit® Neonate (lower body blanket, molex esophageal/rectal probe, set of mittens, set of socks)	5/box
Connecting Hoses (Must purchase one, two or three hoses depending on the kit that is purchased)		
286	9' Reusable Connecting Hose (2.7m)	ea.
286-18	18' Reusable Connecting Hose (5.5m)	ea.
286-27	27' Reusable Connecting Hose (8.2m)	ea.
Hyper-Hypothermia units		
233	Blanketrol® III	115v or 230v

Product Specifications:

Head Wrap Dimensions:

16.4in x 24.6in (41.7cm x 62.5cm)

Patient Vest Dimensions:

53in x 25in (134.62cm x 63.5cm)

Lower Body Blanket Dimensions:

33in x 25in (83.8cm x 63.5cm)

Full Body Blanket Dimensions:

64in x 25in (162.5cm x 63.5cm)
(catalog no. 930 only)



Technology to the next degree™

www.gentherm.com

CLINICAL 24/7 SUPPORT: 1.513.460.2038
TECHNICAL 24/7 SUPPORT: 1.888.437.5608